

Key Feature of the FR3 AED

The key features of this AED are:

- Infant/Child key: converts to a 50 joule setting from 150 joule for adults

- One set of pads for adults and children

- Pads are NOT in a pouch; faster to apply to patient; tested daily.

- SMART CPR – Evaluates key characteristics of the presenting VF and determines the initial therapy: shock first, or CPR first, quickly followed by a shock

- Bilingual configurable (primary language and secondary language)

- Option is to have the AED turn on automatically when the lid is lifted

- Smallest and lightest on the market (good for the bike medics) 3.5 lbs. with battery installed, 8.7" x 5.3" x 2.7"

- Fastest to deliver shock (Less than 8 seconds from end of "stop CPR" prompt.

- Bright, high resolution color screen (2.7" x 2.1")

- Very rugged, crush proof; IP55 rating (jetting water)

- Blue tooth transmission or data card

- 300 shocks or 12 hours of monitoring

- Able to assess the depth of compression and chest recoil; provides coaching to responder (Q CPR)

- Has manual mode so that other rhythms can be shocked if warranted by medical professionals.

- Coaching metronome configurable (as standards change, so can the metronome)

- Training battery converts live AED into a trainer.

- AED Trainer can use expired pads for training

- Daily, weekly, monthly testing with a full joule discharge on battery insertion and monthly tests.